

## CD138 / Syndecan-1 (SDC1) (Plasma Cell Marker) Antibody

Mouse Monoclonal Antibody [Clone SDC1/7177]

| Catalog No      | Format  | Size   |
|-----------------|---|--------|
| 6382-MSM7-P0    | Purified Ab with BSA and Azide at 200ug/ml    | 20 ug  |
| 6382-MSM7-P1    | Purified Ab with BSA and Azide at 200ug/ml    | 100 ug |
| 6382-MSM7-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |


| Applications               | Tested Dillution | Note  |
|----------------------------|------------------|---|
| Immunohistochemistry (IHC) | 1-2ug/ml         | 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes |

### Product Details


|                               |  |
|-------------------------------|--|
| <b>Clone</b>                  | SDC1/7177  |
| <b>Gene Name</b>              | SDC1   |
| <b>Immunogen</b>              | Recombinant fragment (around aa1-250) of human CD138 protein (exact sequence is proprietary) |
| <b>Host</b>                   | Mouse  |
| <b>Clonality</b>              | Monoclonal   |
| <b>Isotype / Light Chain</b>  | IgG1 / Kappa   |
| <b>Mol. Weight of Antigen</b> | 85kDa  |
| <b>Cellular Localization</b>  | Cell surface.  |
| <b>Species Reactivity</b>     | Human  |
| <b>Positive Control</b>       | Human tonsil or lymph node.  |

\*Optimal dilution for a specific application should be determined.


### Product Images for CD138 / Syndecan-1 (SDC1) (Plasma Cell Marker) Antibody



Formalin-fixed, paraffin-embedded human tonsil stained with CD138 Mouse Monoclonal Antibody (SDC1/7177). HIER: Tris/EDTA, pH9.0, 45min. 2°C: HRP-polymer, 30min. DAB, 5min.



Formalin-fixed, paraffin-embedded human tonsil stained with CD138 Mouse Monoclonal Antibody (SDC1/7177). HIER: Tris/EDTA, pH9.0, 45min. 2°C: HRP-polymer, 30min. DAB, 5min.



SDS-PAGE Analysis of Purified CD138 Mouse Monoclonal Antibody (SDC1/7177).  
Confirmation of Purity and Integrity of Antibody.

### Specificity & Comments

CD138, also designated syndecan-1, is a member of the syndecan family of four transmembrane spanning proteins capable of binding to heparan sulfate and chondroitin sulfate molecules. The syndecans main functions are to control cell growth and differentiation as well as to maintain cell adhesion and cell migration. Under normal conditions CD138 is predominantly expressed on mature plasma cells and early preB-cells, while other haematolymphoid cells are negative. Various types of epithelial cells are also CD138 positive. Squamous epithelial cells show strong membranous and some cytoplasmic staining, except for the mature superficial squamous epithelial cells which are unstained. Mature mesenchymal and neural tissues are not stained. Among haematolymphoid neoplasms, CD138 is expressed in practically all cases of plasma cell malignancies. Among non-haematolymphoid neoplasms the expression of CD138 is found in various types of carcinoma.

### Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.

### Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

### Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

### Research Areas

Immunology, Articular Cartilage Extracellular Matrix, B Cell Markers, Cytokine Signaling, Hematopoietic Stem Cells, Infectious Disease, Mesenchymal Stem Cell Differentiation