

NME2 / nm23-H2 / NDPK-B (Suppressor of Metastasis) Antibody

Mouse Monoclonal Antibody [Clone NME2/12061]

| Catalog No | Format | Size |
|------------------|---|--------|
| 4831-MSM10-P0 | Purified Ab with BSA and Azide at 200ug/ml | 20 ug |
| 4831-MSM10-P1 | Purified Ab with BSA and Azide at 200ug/ml | 100 ug |
| 4831-MSM10-P1ABX | Purified Ab WITHOUT BSA and Azide at 1.0mg/ml | 100 ug |


| Applications | Tested Dillution | Note |
|----------------------------|------------------|---|
| Immunohistochemistry (IHC) | 1-2ug/ml | 30 min at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes |

Product Details


| | |
|------------------------|--|
| Clone | NME2/12061 |
| Gene Name | NME2 |
| Immunogen | Recombinant full-length human NME2 protein |
| Host | Mouse |
| Clonality | Monoclonal |
| Isotype / Light Chain | IgG / Kappa |
| Mol. Weight of Antigen | 17kDa |
| Cellular Localization | Cytoplasm, Nucleus |
| Species Reactivity | Human |
| Positive Control | Ubiquitously expressed in all tissues. |

*Optimal dilution for a specific application should be determined.

Product Images for NME2 / nm23-H2 / NDPK-B (Suppressor of Metastasis) Antibody



SDS-PAGE Analysis of Purified NME2 Mouse Monoclonal Antibody (NME2/12061). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human kidney stained with NME2 Mouse Monoclonal Antibody (NME2/12061). Inset: PBS instead of primary antibody; secondary only negative control.

Specificity & Comments

The nm23 gene, a potential suppressor of metastasis, was originally identified by differential hybridization between two murine melanoma sub-lines, one with a high and the second with a low metastatic capacity. Highly metastatic sub-lines exhibit much lower levels of nm23 than less metastatic cells. Based on sequence analysis, nm23 appears highly related to nucleotide diphosphate kinases (NDP). In humans, NDP kinases A and B are identical to two isoforms of human nm23 homologs, namely nm23-H1 and H2, respectively. nm23-H2 is identical in sequence to PuF, a transcription factor that binds to nuclease hypersensitive elements at positions 142 to 115 of the human c-Myc promoter.

Supplied As

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% BSA & 0.05% azide. Also available WITHOUT BSA & azide at 1.0mg/ml.

Storage and Stability

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Research Areas

Immunology

Limitations and Warranty

This antibody is available for research use only and is not approved for use in diagnosis. There are no warranties, expressed or implied, which extend beyond this description. Company is not liable for any personal injury or economic loss resulting from this product.
